TermalVap

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TermalVap is a thermal inhaler equipped with a steel boiler inside of which a 230V electric heating element operates. The appliance generates hot steam which, by way of a sprayer that takes advantage of the Venturi principle, is capable of mixing the steam developed by its own boiler with the main thermal waters located outside of the appliance.

Built with a high thermal insulation plastic body in conformity with European safety standards, the appliance was designed for the inhalation of bottled thermal water and is equipped with specific safety systems such as the double safety thermostat applied on the boiler and the safety cap which is triggered when pressures greater than 2 bar are reached.

The appliance in question should be used as an inhaler for household use.

GENERAL WARNINGS

CAREFULLY READ THE INSTRUCTION MANUAL BEFORE USING THE APPLIANCE. NEVER DISASSEMBLE THE UNIT.

FOR ANY SERVICE CONTACT THE PIC TECHNICAL SERVICE CENTRE. NEVER ALLOW CHILDREN AND / OR INCOMPETENT PEOPLE TO USE THE MEDICAL APPLIANCE WITHOUT PROPER SUPERVISION BY AN ADULT IN FULL POSSESSION OF HIS OR HER MENTAL FACULTIES.

FUNDAMENTAL SAFETY STANDARDS

1. On opening the package ensure that the appliance is whole, taking particular care that there is no damage to the plastic parts which may provide access to the live internal parts of the appliance and ensure that there is no breakage and/or peeling of the power cable. In these cases do not connect the electrical plug and contact the PIC TECHNICAL SERVICE CENTRE.
2. Before connecting the unit always ensure that the electrical data indicated on the silver plate (positioned on the bottom) and the type of plug used correspond with those of the electrical mains you intend to connect to.
3. Before connecting the unit fill the boiler with 800ml of demineralised or distilled water and ensure that there are no water leaks from the bottom of the unit. Otherwise, do not use the appliance, and contact the PIC TECHNICAL SERVICE CENTRE.
4. If there is a visible steam leak from the safety cap (which is an indicator of a malfunction) switch off the appliance and contact the PIC TECHNICAL SERVICE CENTRE. If the appliance does not generate steam after 10 minutes of being switched on there has probably been a drop in power from the mains to which the ap-

pliance was connected. In this case disconnect the appliance from the power socket and carry out the checks as indicated in the MAINTENANCE section. If the problem persists, contact the PiC TECHNICAL SERVICE CENTRE.

5. Observe the indicated safety rules for electrical devices, specifically:
 - Only use original accessories and parts supplied by the manufacturer CA-MI in order to guarantee maximum efficiency and safety of the appliance.
 - NEVER use tap water because the lime scale could compromise proper operation of the appliance, clogging the nebuliser nozzle and the boiler-
 - Check for the presence of demineralised water inside the boiler at every use. Operating the appliance without water will trigger a breaker device which must be manually reset.
 - During use the appliance must be placed on a horizontal and stable surface and the air vents must never be obstructed.
 - Never use the appliance in environments where flammable gaseous mixtures are present.
 - Never submerge the appliance in water.
 - Never touch the appliance with wet and / or damp hands.
 - Do not leave the appliance plugged into the power socket when not in use.
 - Do not pull or tug at the power cable to unplug it from the mains but rather grip the power plug to remove it from the socket.
 - Do not wind the power cord around itself during operation in order to prevent dangerous overheating;
 - Store and use the appliance in environments protected from atmospheric agents and far from any heat sources.
6. This appliance must be used exclusively for the intended use for which it was designed, in other words as a Thermal Inhaler for domestic use. Any use other than the one for which the unit is intended is considered improper and therefore dangerous; The manufacturer will not be liable for any damage caused by incorrect and/or unreasonable use or if the appliance is used in electrical systems that do not conform to prevailing safety regulations. Do not use for purposes other than those foreseen by the manufacturer.
7. Disposal of the medical appliance accessories must be carried out in accordance with the specific prevailing legislation in each country.
8. Keep this manual for further consultation.
9. Disconnect the appliance from the mains, unplugging the plug and switching it off before cleaning and / or maintenance operations.
10. On request the manufacturer will provide electrical diagrams, parts list, descriptions, calibration instructions and / or any other information which may aid technical service personnel in repairing parts of the medical appliance.

11. If the appliance accidentally falls into a sink, bathtub, etc., never try to remove the appliance from the water with the plug inserted in the power socket. Switch off the mains switch, remove the plug from the power socket and contact qualified personnel: do not attempt to operate the appliance until it has been carefully checked by qualified personnel and / or the PIC TECHNICAL SERVICE CENTRE.
12. The medical appliance requires special precautions regarding electromagnetic compatibility and must be installed and put into service following the information provided in the user manual.
13. Some components of the appliance are so small that they can be swallowed by children; therefore keep the appliance out of the reach of children.
14. Keep the accessories out of the reach of children. Children and non self-sufficient people must always use the medical appliance under the close supervision of an adult in full possession of his or her mental faculties. Keep the nebuliser cup out of the reach of children under the age of 36 months in that it contains small parts which can be swallowed. Do not leave the device unattended within the reach of children and/or persons who are not in full possession of their mental faculties as they could strangle themselves with the underwater tube.
15. Keep the appliance well out of the reach of children under the age of 36 months as they could, if unsupervised, strangle themselves with the power cable and / or the underwater tube.
16. In the event of a fault and / or poor operation of the appliance, switch it off and do not tamper with it. **WARNING:** Do not modify this appliance without authorisation from the manufacturer. No electrical and / or mechanical part in the appliance was designed to be repaired by the user. Failure to observe this may compromise the safety of the appliance.
17. **WARNING:** Always position the end of the nebuliser cup where the steam comes out at a distance of 20-40cm from the user's face.
18. Use of the appliance in environmental conditions other than those indicated in this manual may seriously prejudice its safety and technical parameters.
19. In the event of a fault and / or poor operation of the appliance, switch it off and do not tamper with it. For any repairs contact the PIC TECHNICAL SERVICE CENTRE and request the use of original spare parts. Failure to observe this may compromise the safety of the appliance.

Under some fault conditions the housing temperature may be high and may constitute a risk of scalding/burning if these parts are touched. In any case, these temperatures never exceed the limit of 105°C (ref. Interpretation Sheet IEC 60601-1).

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

PIKDARE S.p.A. and CA-MI S.r.L. will not be held liable for accidental or indirect damage should unauthorised modifications, repairs and / or technical operations be carried out on the appliance or if any of its parts have been damaged due to accident, use and / or improper use.

Any operation on the appliance, even if minimum, will immediately invalidate the warranty and in any case does not guarantee correspondence to the technical and safety requirements set forth in directive MDD 93/42/EEC (as updated and amended) and the relative reference standards.

TECHNICAL FEATURES

TYPE (Directive 93/42/EEC)	Class IIa Medical Appliance
MODEL	TermalVap (Thermal Inhaler)
CODE	REF 02 038069 100 000
SAFETY CLASS	Class II
POWER SUPPLY	230 V ~ / 50Hz
POWER ABSORBED	600 W
FUSE	F 1 x 4A L 250V
BOILER CAPACITY	800cc of Distilled water
OPERATION	CONTINUOUS
INHALATION TIME	12 to 15 minutes (for 1,000ml of thermal water)
DIMENSIONS	180 x 250 x 280mm
WEIGHT	1.5Kg (empty)
OPERATING CONDITIONS	Ambient temperature: 5 to 40°C Ambient humidity percentage: 0 to 93% RH without condensation Atmospheric pressure: 700 to 1,060 hPa Altitude: 0 to 2,000 m a.s.l.
STORAGE AND TRANSPORTATION CONDITIONS	Ambient temperature: -25 to 70°C Ambient humidity percentage: 10 to 93% RH without condensation Atmospheric pressure: 500 to 1,060 hPa

SYMBOLS USED ON THE APPLIANCE

	Mark of conformity to directive 93/42/EEC as updated and amended
	General and / or specific warnings
	See the instructions for use
	Class II appliance
	Manufacturer: CA-MI Srl - Via Ugo La Malfa nr.13 Frazione Pilastro - 43013 Langhirano (PR) Italy
	Store in a cool, dry place.
	Storage temperature: -25 to 70°C
	Fuse
	Alternating Current
Hz	Frequency
	Production Batch
	Serial Number

REF	Product identification code	
	On	
○	Off	
IP21	Electrical appliance grade of protection against accidental or intentional contact with the human body and with objects and protection against contact with water.	
	1st DIGIT PENETRATION OF SOLIDS	2nd DIGIT PENETRATION OF LIQUIDS
	Protected against solid bodies larger than Ø 12mm	Protected against vertical falling drops of water
	Atmospheric transportation and storage pressure: 500 to 1,060 hPa	
	Filling the boiler	

A copy of the EC declaration of conformity and a copy of the EC certificate issued by the Notified Body can be requested from CA-MI Srl - Via Ugo La Malfa nr.13 - Frazione Pilastro - 43013 Langhirano (PR) Italy, or alternatively from PIKDARE S.p.A..

The technical specifications may vary without advance notice!

RISKS OF ELECTROMAGNETIC INTERFERENCE AND POSSIBLE REMEDIES

This section contains information on product compliance with the EN 60601-1-2 standard.

The TermalVap model Thermal Inhaler is an electrical medical appliance that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned according to the electromagnetic information provided. Mobile and portable RF communications equipment (mobile phones, transceivers, etc.) may affect the medical appliance and should not be used nearby, adjacent to or on top of the medical appliance. The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as replacement parts, may result in an increase of the appliance or system's emissions or a decrease of its immunity.

The following tables provide information on the EMC (Electromagnetic Compatibility) of this electromedical appliance.

Guide and manufacturer's declaration – Electromagnetic Emissions

The Thermal Inhaler TermalVap ^{REF} 02 038069 100 000 is intended for use in the electromagnetic environment specified below. The customer or user of the TermalVap ^{REF} 02 038069 100 000 must ensure that the device is used in this environment.

Emissions test	Complies	Electromagnetic Environment - Guidance
Radiated/Conducted Emissions CISPR11	Group 1	The TermalVap inhaler ^{REF} 02 038069 100 000 only uses RF energy for its internal function. Therefore its RF emissions are very low and do not cause any interference in nearby electronic equipment.
Radiated/Conducted Emissions CISPR11	Class [B]	The TermalVap inhaler ^{REF} 02 038069 100 000 is suitable for use in all environments, including domestic ones and those directly connected to the public power supply network that supplies buildings used for domestic purposes.
EN 61000-3-2 harmonics	Class [A]	
Voltage fluctuation / flicker EN 61000-3-3	Compliant	

Manufacturer's guide and declaration – Electromagnetic Immunity

The Thermal Inhaler TermalVap ^{REF} 02 038069 100 000 is intended for use in the electromagnetic environment specified below. The customer or user of the TermalVap ^{REF} 02 038069 100 000 must ensure that the device is used in this environment.

Immunity test	Test level	Level of compliance	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) EN 61000-4-2	+/-6kV contact +/-8kV air	The appliance does not change its state	Floors should be wood, concrete or ceramic. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN 61000-4-4	+/-2kV for power supply +/-1kV for signal conductors	The appliance does not change its state	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	+/-1 kV differential mode +/-2 kV common mode	The appliance does not change its state	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations EN 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40 % U_T (60% dip in U_T) for 5 cycles 70 % U_T (30% dip in U_T) for 25 cycles <5 % U_T (>95% dip in U_T) for 5 s	- - -	Mains power quality should be that of a typical commercial or hospital environment. If the TermalVap thermal inhaler user requires the unit to be in continued operation, we recommend powering it from an uninterruptible power supply.
Magnetic Field at mains frequency (50/60 Hz) EN 61000-4-8	3A/m	The appliance does not change its state	The mains frequency magnetic fields must have levels characteristics of installations in commercial or hospital environments.

Note: U_i is the mains voltage

Manufacturer's guide and declaration – Electromagnetic Immunity

The Thermal Inhaler TermalVap ^[REF] 02 038069 100 000 is intended for use in the electromagnetic environment specified below. The customer or user of the TermalVap ^[REF] 02 038069 100 000 must ensure that the device is used in this environment.

Immunity test	Level indicated by IEC 60601-1-2	Level of compliance	Electromagnetic Environment - Guidance
Conducted immunity EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	$V_1 = 3 \text{ V rms}$	Portable and mobile RF communication devices must not be used closer to any part of the TermalVap device, including the cables, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance
Radiated immunity EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	$E_1 = 3 \text{ V / m}$	$d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P} \text{ from 80 MHz to 800MHz}$ $d = [23 / E_1] \sqrt{P} \text{ from 800 MHz to 2.5GHz}$ <p>Where P is the maximum rated output power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). The intensities of the fields from fixed RF transmitters, as determined in an electromagnetic study of the site^{a)}, may be less than the level of conformity of each frequency interval^{b)}. Interference may occur near devices marked by the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz the range of the highest frequency is applied.
 Note 2: These guidelines may not cover all situations. Electromagnetic radiation is affected by absorption and reflection by structures, objects and people.
 a) The field intensities for fixed transmitters like base stations for radio telephones (cellular and cordless) and terrestrial mobile radios, HAM radio devices, AM and FM radio transmitters and TV transmitters can not be theoretically foreseen accurately. In order to establish an electromagnetic environment caused by fixed RF transmitters an electromagnetic study of the site must be considered. If the field intensity measured in the location where the appliance is used exceeds the applicable level of conformity above the normal operation of the appliance itself must be placed under observation.
 If abnormal performance is noted additional measures may required such as different orientation or positioning of the appliance.
 b) The field intensity on a frequency range from 150 kHz to 80 MHz must be less than 3 V/m.

Recommended separation distances between portable and mobile radio communication devices and the electromedical appliance

The Thermal Inhaler TermalVap ^[REF] 02 038069 100 000 is intended for use in an electromagnetic environment where RF radiated disturbances are under control. The TermalVap ^[REF] 02 038069 100 000 customer or operator may contribute to preventing electromagnetic interferences by ensuring a minimum distance between mobile and portable RF communication devices (transmitters) and the TermalVap ^[REF] 02 038069 100 000 appliance and recommended below in relation to the maximum output power of the radio communication devices.

Maximum rated output power of the transmitter W	Separation distance at the transmitter frequency m		
	150 kHz at 80 MHz $d = [3.5\sqrt{V_i}] \sqrt{P}$	80 MHz at 800 MHz $d = [12/E_i] \sqrt{P}$	800 MHz at 2.5 GHz $d = [23/E_i] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters with maximum rated output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the frequency of the transmitter where P is the maximum rated output power of the transmitter in Watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the distance for the highest frequency range is to be applied.

Note 2: These guidelines may not cover all situations. Electromagnetic radiation is affected by absorption and reflection by structures, objects and people.

CA-MI S.r.l. declares to have verified the reciprocal compatibility of the components which make up the medical appliance and to have completed the assembly in accordance with the instructions contained in the instruction manual.

The adjustable nebuliser cup (1) and the underwater tube (8) must be replaced if there are visible signs of deterioration of the material that makes up these accessories.

Life expectancy: 200 treatment cycles (1 treatment = 12 applications).

Water tank cap: we recommend replacing the valve cap every 3 years or every 600h of operation.

CLEANING OPERATIONS OF THE MAIN UNIT

To clean the external part of the appliance use a damp cotton cloth with detergent. Do not use abrasive detergent substances and solvents.

TAKE PARTICULAR CARE TO ENSURE THAT THE INTERNAL LIVE PARTS OF THE APPLIANCE DO NOT COME INTO CONTACT WITH LIQUIDS. NEVER WASH THE APPLIANCE UNDER WATER OR BY SUBMERGING IT.

Disconnect the appliance from the mains, unplugging the plug and switching it off before cleaning and / or maintenance operations.

CLEANING THE ACCESSORIES

Carry out the cleaning operations as indicated below.

Cleaning the steam conveying nebuliser cup (1):

Extract the steam conveying nebuliser cup from the nebuliser cup nozzle seat (2) unscrewing the fixing screw (4). Submerge the nebuliser cup in a common anti-lime scale liquid and allow it to work until the lime scale deposits disappear completely. Observe the anti-lime scale product manufacturer's instructions. Rinse generously under running water.

Cleaning the sprayer nozzle:

Extract the steam conveying nebuliser cup (1) from the nebuliser cup nozzle seat (2) unscrewing the fixing screw (4). Remove the nozzle guard and, using the included CH17 wrench, unscrew the nut that secures the component to the boiler. Use your hand to grip the plastic body of the nozzle and unscrew in an anticlockwise direction. To prevent lime scale encrustations from forming fill a common cup with 50% water and 50% vinegar. Submerge the sprayer in the prepared solution for 1 hour and then rinse under running water. To remove lime scale encrustations from the sprayer fill a common cup with 50% water and 50% common anti-lime scale solution. Submerge the component in the prepared solution for 12 hours and then rinse very well under running water. In order to remove encrustations from the centre hole of the sprayer you can also use the included pin (14) periodically. If the sprayer was disassembled for cleaning, correctly reassemble it: the nut gasket must be up against the bush (5) that comes out alongside the carcass.

BEFORE AND AFTER EACH USE AND CLEANING, ASSEMBLY AND DISASSEMBLY OPERATION THE INTEGRITY OF EACH PART OF THE PRODUCT AND THE ACCESSORIES USED MUST BE CHECKED. CAREFULLY CHECK THAT THE SPRAYER WAS CORRECTLY SECURED TO THE BODY OF THE APPLIANCE AND IN PARTICULAR THAT THE PLASTIC NEBULISER CUP IS WELL SECURED AND WHOLE IN ALL ITS PARTS (NO CRACKS OR BREAKS, INTERNAL PROTECTION WALL PRESENT). BREAKS OR INCORRECT ASSEMBLY MAY GENERATE DANGEROUS CONDITIONS DURING NORMAL AND CORRECT OPERATION.

INSTRUCTIONS FOR USE

1. The appliance should be checked before every use in order to detect any operational faults and / or damage due to transport or storage.

2. During inhalation the patient must sit in an upright and relaxed position at a table and not in an armchair, positioning the appliance on a flat and stable surface.
3. Remember not to hold the appliance in your hand during treatment and / or to avoid prolonged body contact with the appliance.
4. Unscrew the safety valve cap (6) (the rotation direction is indicated with the letter A) and fill the boiler with 800ml of demineralised or distilled water using the graduated decanter (10) and funnel (13) which are included. The boiler filling operations must be carried out with the appliance not powered. Screw the valve cap back in completely, tightening it securely (the rotation direction is indicated with the letter C).

During normal operation do not remove the guard from the sprayer. NEVER submerge the appliance under a direct stream of running water. For the boiler filling operations always use the graduated condensation collection decanter and the funnel which are included.

5. Insert the nebuliser cup (1) in the sprayer and turn the stop (4) to lock it in place. The nebuliser cup depleted water drain must be positioned downward.
6. Place the decanter (10) under the nebuliser cup drain (3) in order to collect the discharged thermal water.
7. Take a bottle of thermal water and position it in the specific compartment. Then position the underwater tube (8) in the same bottle. Ensure that the tube is positioned on the bottom of the bottle and close the bottle with the closing cap (12).
8. At the other end of the underwater tube connect it to the sprayer connector (11).
9. Before proceeding with activation of the appliance ensure that you have filled the boiler with 800 ml of distilled or demineralised water and that you have completely screwed the valve cap (6) back on.
10. Connect the appliance electrical plug to the electrical mains.
11. To begin treatment switch on the appliance by pressing the luminous switch (7) to the I position. The green switch will light up to indicate that the appliance is correctly powered. Wait for an abundant steam jet emission (about 6 – 8 minutes), indicating that the device is ready for use.
12. Sit in front of the appliance at a distance of 20-40 cm and begin the treatment which will last about 12-15 minutes (for 1L of thermal water).

TAKE PARTICULAR CARE NOT TO REMOVE THE BOILER GUARDS PLACED NEAR THE NEBULISER NOZZLE

13. Once the treatment is over position the switch to 0 and remove the plug from the mains socket; Before opening the valve and emptying the boiler wait until the appliance is completely cooled. Then, once the bottle has been removed, remove the accessories (separate the sprayer nebuliser cup and the underwater tube), unscrew the safety cap and empty the residual water from the boiler.
14. **WARNING:** The quantity of water introduced into the boiler (800cc) allows only one inhalation to be carried out. The quantity of depleted water (collected in the decanter) has no more therapeutic value and must therefore not be reused.
15. Proceed with the cleaning operations as indicated in the CLEANING section.
16. Check the water level inside the boiler at every treatment.

During operation of the device high temperatures are registered on the steam conveying nebuliser cup because it is close to the aqueous steam emission point needed for treatment and for mixture with the thermal water. In any case this temperature does not constitute a risk as the treatment must be carried out maintaining a distance of about 20-40cm: therefore, at this distance these temperatures are not dangerous.

Never come into contact with the steam conveying nebuliser cup during operation of the device.

DO NOT OPERATE THE APPLIANCE WITHOUT WATER IN THE BOILER AS THIS WILL TRIGGER A THERMAL BREAKER THAT MUST BE RESET MANUALLY.

THIS THERMOSTAT MUST BE RESET AT AN AUTHORISED CENTRE AND / OR THE PIC TECHNICAL SERVICE CENTRE.

WARNING: The power cable is the separating element from the electrical mains; even though the appliance has a specific on / off switch the power plug must be kept accessible when the appliance is in use in order to allow a further method if necessary for disconnection from the mains.

THERMAL BREAKER: The appliance has a thermal breaker that trips if excessively high temperatures are reached. If this occurs, switch off the appliance by pressing the luminous switch (toward the 0 position) and allow it to cool for about 30 minutes. The thermal breaker may trip for one of the following reasons:

- the appliance operated outside of the operating limits indicated in this manual.

- the appliance operated close to heat sources and / or in environments with high temperatures.

The appliance itself also has a second thermal safety device which must be manually reset that trip in the event the first one is faulty, definitively locking up the appliance. The second thermostat also trips if no water is placed in the boiler, preventing the occurrence of high temperatures. In this case, in order to reuse the appliance contact the PIC TECHNICAL SERVICE CENTRE which will reset the safety thermostat after checking that the system still works.

PERIODIC MAINTENANCE CHECK

The TermalVap appliance has no parts that require maintenance and/or lubrication (except for cleaning operations). Nevertheless, some simple checks must be carried out in order to verify operation and safety of the appliance before each use. Remove the appliance from the box and always check to ensure that the housing and the power cable are whole as they could have been damaged during the previous use. After filling the boiler with water, screw on the valve cap and ensure that the rubber gasket is up against the boiler bush which comes out of the appliance's carcass. If the sprayer was disassembled for cleaning, ensure that it has been correctly reassembled. The nut gasket must be up against the bush that comes out alongside the carcass. Then connect the cable to the mains and switch on the switch; Ensure that there are no steam leaks from the valve cap and/or the sprayer nut. At the end of inhalation always cool and remove the water from the boiler (see the INSTRUCTIONS FOR USE section). The appliance is protected by a fuse (1 x F4A L 250V) located inside it which cannot be replaced from the outside, therefore you must contact qualified personnel to have it replaced.

Defect type	Cause	Solution
1. Water is not drawn from the bottle	Sprayer clogged with lime scale	Let the appliance cool, disassemble the sprayer with the included wrench and clean as indicated in the cleaning section.
2. The appliance switches off at mid inhalation	Sprayer clogged with lime scale	Let the appliance cool, disassemble the sprayer with the included wrench and clean as indicated in the cleaning section.
3. Steam is not produced	The power supply cable is not connected. The safety thermostat has tripped	Ensure that the appliance's power cord is properly connected to the electrical mains and that the switch is set to the I position Ensure that there is water in the boiler. If the safety thermostat has tripped contact the PIC TECHNICAL SERVICE CENTRE.
4. Steam comes out at low pressure from the nebuliser nozzle	Nebuliser nozzle blocked	Let the appliance cool, disassemble the sprayer with the included wrench and clean as indicated in the cleaning section.
5. The appliance switches off without completely nebulising the water in the bottle	Nebuliser nozzle blocked	The nebuliser nozzle may be blocked. For cleaning see the "Cleaning Accessories" section
6. The appliance switches off during inhalation	Safety thermostat tripped	The safety thermostat may have tripped. Wait 30 for the appliance to cool. If the problem persists contact the PIC TECHNICAL SERVICE CENTRE.
Defects 1 - 2 - 3 - 4 - 5 - 6	None of the solutions worked	Contact the PIC TECHNICAL SERVICE CENTRE

In the event of a fault, contact the PIC SERVICE CENTRE and / or authorised personnel. Do not open the appliance under any circumstances. If technical support personnel must replace the power cable always ensure correct connection and adequate anchoring of this component.

BEFORE CARRYING OUT ANY CHECKS, IN THE EVENT OF ABNORMALITIES OR MALFUNCTIONS CONTACT THE PIC TECHNICAL SERVICE CENTRE. PIKDARE S.P.A. AND CA-MI S.R.L. DO NOT OFFER ANY TYPE OF WARRANT

TY FOR APPLIANCES WHICH ARE FOUND TO HAVE BEEN TAMPERED WITH UPON INSPECTION BY TECHNICAL SERVICE.

CONDITIONS OF VALID WARRANTY - REPAIR METHODS

IN OBSERVANCE OF THE NEW EUROPEAN REGULATIONS, PIKDARE S.P.A. AND CA-MI S.R.L. LIST A FEW FUNDAMENTAL POINTS FOR PRESERVING THE HYGIENE OF THE APPLIANCE AND THE OPERATORS WHO USE IT. CA-MI S.R.L. AND PIKDARE S.p.A. TRUST THAT THESE REGULATIONS WILL BE OBSERVED IN ORDER TO GUARANTEE THE HYGIENE AND HEALTH OF ALL OPERATORS IN ORDER TO OBTAIN QUALITY AND WELLNESS.

Any appliance which is returned to the PIC Technical Service Centre will be subjected to hygiene checks before being repaired. If the PIC Technical Service Centre deems the appliance unsuitable for repair due to visible signs of external and/or internal contamination, the appliance will be returned to the customer clearly marked APPLIANCE NOT REPAIRED with a letter of explanation on the defects found attached to it. CA-MI S.r.l. and the PIC Technical Service Centre will assess whether the contamination is caused by poor operation or incorrect use. CA-MI S.r.l. and PIKDARE S.p.A. will not be liable for any accessories that show signs of contamination, therefore the PIC Technical Service Centre will replace them and charge the material cost to the customer.

As concerns the above, it is COMPULSORY to carefully disinfect the external carcass using a cloth dampened with denatured alcohol or hypochlorite based solutions, as well as the accessories, submerging them in the same disinfecting solutions. Place them in a bag labelled disinfected appliance and accessories. We ask that you always specify the defect found in order to be able to carry out the repair as quickly as possible. We therefore ask that you carefully read the instructions for use in order to prevent compromising the appliance with unsuitable use. We ask that you always specify the defect found so that the technicians at the PIC Technical Service Centre will be able to assess whether or not the defect is covered under warranty. The warranty is valid for 24 months from the date of purchase against any material or construction defect as long as no unauthorised operations or tampering has been carried out on the appliance. The warranty is valid only if the appliance is returned to the PIC TECHNICAL SERVICE CENTRE accompanied by the fiscal receipt which proves the purchase date. The term warranty means the replacement or repair of the appliance's components which are found to be defective from their origin due to manufacturing defects.

The Warranty does not cover:

- The included accessories and parts subject to normal wear.
- Operations on presumed defects or for convenience.
- Operations for break downs or defects due to improper use, neglect or negligence and if the damage is not imputable to the manufacturer (accidental dropping, careless transportation, etc.).

The warranty does not imply any reimbursement for direct or indirect damages of any nature to persons, animals or objects stemming from the improper use of the product or its period of inefficiency.

The warranty is valid from the product date of purchase as certified by the fiscal receipt.

CA-MI Srl
Via Ugo La Malfa nr.13 - Frazione Pilaastro
43013 Langhirano (PR) Italia
Tel (0039) 0521 637133 - 631138
Fax (0039) 0521 639041
E-mail: vendite@ca-mi.it; export@ca-mi.it
Web Site: www.ca-mi.it

Distribuito da PIKDARE S.p.A.
Via Saldarini Catelli, 10
22070 Casnate con Bernate (CO)-Italy
Made in Italy

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Numero Verde
800-030-911
Linee libere 08:00 - 18:00 e 14:00 - 17:00

www.picsolution.com

Numero verde
consumatori Italia
800 900 080

www.picsolution.com

REV. 01/2019

REF 02 038069 100 000