Termometro auricolare ad infrarossi

Comfort QUICK

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REF 00000656000000

Radiant Innovation Inc.
Indirizzo: 1F, n.3, Industrial East 9th Rd, Science-Based Industrial Park, HsinChu, Taiwan

MT Promedt Consulting GmbH
Indirizzo: Altenhofstrasse 80, D-66386 St. Ingbert, Germania

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EN Infrared ear thermometer

Comfort Quick

Dear Customer, thank you for choosing our product!

This thermometer was designed by a panel of highly proficient technicians and produced in conformity with all applicable European Standards. It uses the smallest probe in the world (5.9 mm in diameter) and is therefore ideal for measuring ear temperatures of newborns, infants and all children in general.

The thermometer is an adjusted mode thermometer that picks up the infrared radiation emitted by the tympanic membrane and converts the reading into an equivalent oral measurement.

Always keep the thermometer and probe caps out of reach of children and never leave children unattended when using the thermometer. The thermometer contains small parts that children can accidentally swallow. If in doubt, or if the temperature goes up, contact your doctor for advice.

Please report to the manufacturer and the competent authority of the Member State in which you are established about any serious incident that has occurred in relation to this device.

BODY TEMPERATURE

It is important to bear in mind that body temperature varies from person to person, depending on where the reading is taken, and it varies throughout the day, fluctuating between 35.5 °C and 37.8 °C (95.5 °F -100 °F). This is why, in order to evaluate the temperature reading properly, we recommend you record the child's normal ear temperature, when he is in good

health, at various times of day and use this as your benchmark when assessing if the temperature has risen or not. To determine whether someone has a fever, it is recommended to compare the current temperature reading with one taken in good health: an increase equal to or greater than 1°C (1.8 °F) compared to the benchmark temperature is usually indicative of a fever. Furthermore, not only do different measurement points (rectal, underarm, oral, frontal, ear) provide different readings, but it is also difficult to produce valid correlations for each subject (there are only statistical data in this sense). It is hence absolutely normal that the temperature differs when taken in the ear rather than the underarm and it is wrong to try and compare the two. It is also wrong to compare readings taken on the right and left sides of the body as the temperature readings depend on the depth of the blood vessels where the reading is taken, and this depth usually differs between the right and left side. The indicative readings below refer to standard adult temperatures taken in different parts of the body:

- ear temperature: 35.8 – 38°C / 96.4 – 100.4° F
- rectal temperature: 36.6 – 38°C / 97.9 – 100.4° F
- underarm temperature: 34.7 – 37.3°C / 94.5 – 99.1° F
- oral temperature: 35.5 – 37.5°C / 95.9 – 99.5° F

HOW TO TAKE A TEMPERATURE

Precautions

- The thermometer tends to compensate ear temperature readings according to the room temperatures. To get the most accurate reading, it is therefore recommended to place the thermometer in the room where it will be used for at least 20-25 minutes. The person whose tem-

perature will be taken should also be in the same room for at least 5 minutes before taking the reading. Whilst taking the temperature, make sure that the patient and thermometer are not exposed to direct heat, sunlight or air conditioning.

- Before taking the reading, check that the ear canal is clean and accessible; do not clean the ear just before using the probe as the use of hot or cold water will affect the reading.
- Do not take a person's temperature for at least thirty minutes after physical effort, a bath or after a meal.
- If the probe is pointed at a hot object, it will perform a remote temperature reading. To avoid such reading errors, keep the time from when it is switched on to when it takes the ear reading to a minimum.
- Do not touch the sensor, the probe or the probe cap with your fingers. If fingerprints get left on the sensor, remove them following the instructions provided in the "CARE AND MAINTENANCE" section.

READING PROCEDURE

Important: before switching the appliance on, the probe cap must be fitted correctly, following the instructions provided in the "Fitting the probe cap" section.

Fitting the probe caps.

1. Place the probe cap in the housing in the thermometer case (see diagram 1). Important: position the probe cap making sure the adhesive part, i.e. the knurled section, (see diagram 1) is facing upwards.
2. Align the centre of the probe cap with the centre of the probe.
3. Insert the probe into its seat until you hear it click into place; this

mean it is now ready to use.

Fig. 1

If the probe cap is not fitted correctly, when the appliance is switched on, the ► symbol will blink on the display and it will not be possible to take any readings (the appliance will beep three times and no message will appear on the display).

Fig. 2

If this happens, fit the probe cap again following the instructions provided in points 1-3 of this section very carefully.

Warning! Only use dedicated probe caps [REF](#) 00016395000000

Switching on the appliance and taking readings

Press the  (ON/MEM) button on the left side of the appliance.

Fig. 3

When the appliance is switched on, it will run an autotest and all the segments on the display will light up.

If any segments do not light up, do not use the thermometer under any circumstances.

When the autotest sessions has been completed, the thermometer display will look like the picture in diagram 4

The thermometer is ready for use as soon as the  (ear) icon stops blinking

and the appliance beeps twice.

Fig. 4

If an error message appears on the display, please refer to page 19.

Important: Before taking a reading, always make sure that the sensor is clean and has not been damaged in any manner; also check that the ear canal is perfectly clean too. Now press the  (ON/MEM) button to turn the appliance ON.

How to take a temperature:

1. Gently move the ear back slightly to straighten the ear canal and very carefully insert the probe into the canal. Aim towards the timpani membrane to get the most accurate reading.

Diag. 5

- For children under two years of age, gently pull the ear backwards.
- For children over two years of age and adults, gently pull the ear upwards and backwards.

0-2 years

>2 years

Fig. 6

2. Press the "SCAN" button on the right side of the thermometer and in about 1 second the thermometer will emit an acoustic signal which means it has taken the reading. The temperature reading will be visible

on the display immediately.

Fig. 7

Note: Ear temperatures have to be converted into the equivalent oral temperature.

The thermometer converts the temperature to show the equivalent oral temperature (depending on the result of clinical Bias).

Note: the reading taken in the right ear canal can vary to the reading in the left ear canal: this is because the depth and size of the right and left tympanic artery can be different.

Once the temperature has been recorded, the thermometer will emit another short beep when it is ready to be used again.

IMPORTANT

Always take temperatures in the same ear canal to obtain more accurate results. Always use the same ear if the temperature has to be taken repeatedly at various intervals. If consecutive readings are required, wait at least 2 minutes between readings to prevent the sensor from overheating and recording inaccurate measurements.

Fig. 8

Holding the thermometer in your hand for too long and/or exposing it to heat sources could cause the sensor to record a higher than normal room temperature. This could cause the ear temperature reading to be lower than normal.

"LAST READING" FUNCTION

The  icon indicates that the "last reading" function is enabled. When

this function is used, you can take temperatures one after another, without switching the thermometer off, and the last reading will be visible in the top right corner of the display.

MEMORY FUNCTION

This thermometer has a memory which can store up to 25 readings (Please note: if repeated readings are taken without switching the thermometer off in between, the memory will only store the last reading taken). To access the memory, switch the thermometer on, press the  (ON/MEM) button repeatedly to scroll through the various readings until you find the one you need.

CHANGING BETWEEN FAHRENHEIT (°F) AND CELSIUS (°C)

To change the reading mode on the LCD display from °C to °F: switch the appliance off, press and hold the "SCAN" button, now press and release the  (ON/MEM) button. When the "°F" sign appears on the LCD display, release the "SCAN" button. Use the same procedure to change the reading mode on the LCD display from °F to °C.

FEVER INDICATION

If the thermometer detects a reading of $\geq 37.5^{\circ}\text{C}$ (or 99.5°F) it will emit an acoustic signal followed by three short beeps to warn the user that the patient could have a fever.

SWITCHING OFF THE THERMOMETER

The appliance switches off automatically if unused for over 1 minute to preserve battery life, or if the  but-

ton is pressed and held for 5 seconds.

CARE AND MAINTENANCE

The sensor (shown in the diagram 9) is the most delicate part of the thermometer.

Fig. 9

Handle with great care when cleaning the sensor to avoid damaging it in any way.

- Store the appliance in a dry place, away from liquids and direct sunlight.
- Use cotton wool or pads soaked in 70% rubbing alcohol to clean the sensor and the plastic part which comes into contact with the ear canal.
- Use a new probe cap for each reading to avoid any cross contamination. Clean the sensor if dirty to ensure accurate readings. Allow the sensor to dry completely for at least 1 minute.

Important notes:

- Check that the sensor is clean.
- before taking someone's temperature, make sure the person is in the room for about 5 minutes before the reading and has not undergone physical exertion or had a bath in the last 30 minutes; place the thermometer in the room at least 20/25 minutes before the reading so it can adapt to the room temperature.
- Use a new probe cap for each reading.
- Use dedicated probe caps only
- Check that the ear canal is perfectly clean before taking someone's temperature. Try to avoid cleaning it just

before the reading

6. Keep a distance from all sources of heat when taking temperatures.
7. The first time you use the thermometer, remember to remove the battery protection tab (see paragraph "HOW TO INSERT/REPLACE THE BATTERY")
8. If the thermometer or the probe get damaged in any manner, the thermometer could encounter difficulties whilst taking temperatures: if the thermometer is dropped or there is sign of damage, do not use it again before having it checked
9. Have the calibration checked every 3 years. For further information please consult the Chicco website at www.chicco.com


⚠ Choking from swallowing probe covers and connection ring by children or pets is possible, please keep the probe covers and connection ring at places where children and pets can't reach.

IMPORTANT: this thermometer is NOT waterproof. Always store the appliance in a clean, dry place and at a distance from liquids. Storage temperature: -20°C ~ 50°C (-4°F ~ 122°F)

HOW TO INSERT/REPLACE THE BATTERY

The appliance comes with a lithium battery (CR2032x1 pc).

Fig. 10

1. Insert a small pointed object in the hole on the back of the thermometer whilst applying pressure on the battery seat which will then slide

out easily (see diagram 10a).

2. When using the appliance for the first time, remember to remove the battery safety tab.
3. To replace a spent battery, hold the thermometer firmly and insert a screwdriver under the battery to remove it (see diagram 10b)
4. Insert the new battery under the metal hook on the left hand side (number 1), as seen in diagram 10c and press the battery downwards (number 2) until you hear it click into place. Important: insert the battery with the positive pole (+) facing upwards.
5. Replace the battery cover.

Warning: replace with a CR 2032 battery only.

Never use other types of batteries to prevent fire and explosion risks.

Warning: the battery can explode if not handled with care. Do not recharge, pull apart, overheat or set the battery on fire.

⚠ **Keep the battery out of reach of children. To dispose of the thermometer, please follow the instructions provided by EC Directive 2012/19/EC (see page 28); to dispose of the battery, please follow the instructions provided by EC Directive 2006/66/EC (see page 28).**

TECHNICAL SPECIFICATIONS

1. Complies with EN 60601-1 and EN 60601-1-2 Standard requirements.
2. Complies with EN 12470-5 Standard requirements for Clinical Thermometers – Section 5: IR ear thermometer performances (with maximum appliance)
3. Temperature reading range:

- 34~42.2°C (93.2~108°F)
4. Accuracy level: $\pm 0.2^{\circ}\text{C}$ (0.4°F) between $35\sim 42^{\circ}\text{C}$ ($95\sim 107.6^{\circ}\text{F}$) and $\pm 0.3^{\circ}\text{C}$ (0.5°F) beyond this range
 5. Operating range: $10\sim 40^{\circ}\text{C}$ (50°F ~ 104°F); relative humidity 15/85 %
 6. The "Repeatability" is 0.17°C .
 7. Battery life: around 3,000 readings. Enable normal use: 1 year.
 8. Expected Service Life: 4 years.
 9. Enclosure Rating: IP22
 10. Dimensions: 120.4 x 37.8 x 40.3 mm; Weight: 57.7 grams including battery.

- A. Storage temperature range: $-20\sim +50^{\circ}\text{C}$, RH $\leq$ 85%
- B. Transportation temperature shall be less than 70°C , RH $\leq$ 95%
- C. Atmospheric pressure: 800~1013 hPa.

 The device should not submerge into any liquids and expose it to direct moisture.

INTENDED USE

The ear thermometer is electronic thermometer using an infrared detector (thermopile detector) to detect body temperature from ear canal in people of all ages.

INTENDED OPERATOR

Have eight years of education, no maximum.

Ex.SN:E512A000001, the first "E" is External, the second number "5" is the last number of manufacture year, the third and the fourth number "12" is the manufacture month, the others is the serial number.

Classified to EN 60601-1 International Standards :

1. Internally powered appliance.
2.  Applied part: BF type.
3. Continuous operation mode

CE 1639: conforms to the 93/42/CEE + 2007/47/EC Medical Equipment Standard.

Note:

1. The thermometer was produced in compliance with all standards currently in force and meets all magnetic compatibility protection requirements. However, it is necessary to take certain EMC precautions and prepare and use the thermometer following the EMC instructions provided in this document.

2. Portable and mobile radio communication devices can affect the correct functioning of electro-medical devices. In the presence of electro-magnetic disturbance, it is highly.

Manufacturer's declaration-electromagnetic emissions		
The Comfort Quick REF 00000656000000 is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the Comfort Quick REF 00000656000000 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance (for home healthcare environment)
RF emissions CISPR 11	Group 1	The Comfort Quick REF 00000656000000 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Comfort Quick REF 00000656000000 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Manufacturer's declaration-electromagnetic immunity			
The Comfort Quick REF 00000656000000 is intended for use in the electromagnetic environment (for home healthcare) specified below. The customer or the user of the Comfort Quick REF 00000656000000 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance (for home healthcare environment)
Electrostatic discharge(ESD) IEC 61000-4-2	Contact:±8 kV Air±2 kV,±4 kV,±8 kV,±15 kV	Contact:±8 kV Air±2 kV,±4 kV,±8 kV,±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Power frequency(50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz and 60 Hz	The Comfort Quick REF 00000656000000 power frequency magnetic fields should be at levels characteristic of a typical location in a typical home healthcare environment.

Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^{a)} should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 - At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^{a)} Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Comfort Quick [REF] 00000656000000 is used exceeds the applicable RF compliance level above, the Comfort Quick [REF] 00000656000000 should be observed to verify normal operation. If abnormal performance is observed, additional measures <u>may</u> be necessary, such as re-orienting or relocating the Comfort Quick [REF] 00000656000000.			

Manufacturer's declaration-electromagnetic immunity Test specifications for ENCLOSURE PORT IMMUNITY to RF wire- less communications equipment							
The Comfort Quick REF 00000656000000 is intended for use in the electroma- gnetic environment (for home healthcare) specified below. The customer or the user of the Comfort Quick REF 00000656000000 should assure that it is used in such an environment.							
Test fre- quency (MHz)	Band ^{a)} (MHz)	Servi- ce ^{a)}	Modula- tion ^{b)}	Maxi- mum power (W)	Distance (m)	IMMUNI- TY TEST LEVEL (V/m)	Complian- ce LEVEL (V/m) (for home healthca- re)
385	380 – 390	TETRA 400	Pulse modula- tion ^{b)} 18 Hz	1,8	0,3	27	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ±5 kHz devia- tion 1 kHz sine	2	0,3	28	28
710	704 – 787	LTE Band 13, 17	Pulse modula- tion ^{b)} 217 Hz	0,2	0,3	9	9
745							
780							
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modula- tion ^{b)} 18 Hz	2	0,3	28	28
870							
930							
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modula- tion ^{b)} 217 Hz	2	0,3	28	28
1 845							
1 970							

2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{a)} 217 Hz	2	0,3	28	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{a)} 217 Hz	0,2	0,3	9	9
5 500							
5 785							

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Recommended separation distance between portable and mobile RF communications equipment and the Comfort Quick ^{REF} 00000656000000.

The Comfort Quick ^{REF} 00000656000000 is intended for use in an electromagnetic environment (for home healthcare) in which radiated RF disturbances are controlled. The customer or the user of the Comfort Quick ^{REF} 00000656000000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Comfort Quick ^{REF} 00000656000000 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 KHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,7 GHz $d = 2,3 \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 - At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Symbol Descriptions					
	The CE mark and Notified Body Registration Numbers, the requirement of Annex II from Medical Device Directive 93/42/EEC are met.		It means that the device is subject to the EU directive on the disposal of electrical and electronic equipment. To protect the environment, dispose of useless device at appropriate collection sites according to national or local regulations.		Authorized representative in the European community
	Caution		Instructions		Paper Recycling
	Manufacturer		BF type applied part		Battery Recycling
	Classification for water ingress and particulate matter.		Date of production		Medical device

WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use.

The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events.

For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

Troubleshooting:

Error message	Problem	Action
	The probe cap is not fitted correctly	Check that the probe cap is positioned correctly (see the "Fitting the probe cap" section).
	Error 5-9, the system is not working properly.	Remove the battery, wait one minute, replace it and switch it back on. If the message appears again, please contact your retailer for further assistance.
	The battery has run low and cannot take any temperature readings.	Replace the battery.
	Readings taken before the appliance is ready.	Wait for the icons to stop blinking before using the appliance
	The temperature reading is not within the foreseen range of from 10°C to 40°C (50°F - 104°F).	Place the thermometer in a room for at least 30 minutes with a temperature of from: 10°C to 40°C (50°F-104°F).
	The temperature reading is over 42.2 °C (108.0°F).	Make sure the patient remains in the same room for about 5 minutes and has not undergone physical exertion or had a bath in the last 30 minutes; then check that the probe is clean before taking the temperature again.
	The temperature reading is under 34°C (93.2°F) in forehead mode.	If it still malfunctions, please contact your local retailer.
	The appliance will not switch on.	Replace the battery with a new one. If it still malfunctions, please contact your local retailer.

If the LCD displays the following icons at intervals:

this means that the appliance is in Calibration mode. The readings are not accurate when the appliance is in this mode. Switch off the appliance and switch it back on again to ensure it functions properly.

REF 00000656000000

Radiant Innovation Inc.
Address: 1F, n.3, Industrial East 9th Rd,
Science-Based Industrial Park,
HsinChu, Taiwan

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Read these instructions carefully and
keep them for future reference.

**This product complies
with the Directive
2012/19/EU.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

**THIS PRODUCT COMPLIES
WITH EC DIRECTIVE
2006/66/EC**

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate

waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Led. The user is responsible for taking the batteries to a special waste disposal site at the end of their life, so that they can be treated and recycled. If the spent batteries are collected correctly as separate waste, they can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and human health, and contributes towards the recycling of the batteries' substances. Non-compliance with the norms on battery disposal damages the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the batteries.

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Artsana S.p.A.
Artsana S.p.A. - Via Saldarini Catelli 1 - 22070 Grandate (CO) Italy